

Product datasheet for **RA19014-100**

GABARAP Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Purified
Gene Name:	Homo sapiens GABA(A) receptor-associated protein (GABARAP)
Database Link:	Entrez Gene 56486 MouseEntrez Gene 58974 RatEntrez Gene 11337 Human
Synonyms:	GABA Receptor Associated Protein, FLC3B
Protein Families:	Druggable Genome
Protein Pathways:	Regulation of autophagy



[View online »](#)