

Product datasheet for RA19044-100

Prostate Specific Antigen (KLK3) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens kallikrein related peptidase 3 (KLK3), transcript variant 3
Database Link:	Entrez Gene 354 Human
Synonyms:	Prostate-specific antigen, Semenogelase, Gamma-seminoprotein, P-30 antigen
Protein Families:	Druggable Genome, Protease, Secreted Protein
Protein Pathways:	Pathways in cancer, Prostate cancer



[View online »](#)