

Product datasheet for **RA19057-100**

Shprh Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Mus musculus SNF2 histone linker PHD RING helicase (Shprh), transcript variant 1
Database Link:	Entrez Gene 268281 Mouse
Synonyms:	KIAA2023



[View online »](#)