


Product datasheet for RA19065-100


Grm1 Rabbit Polyclonal Antibody

Product data:



Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Rattus norvegicus glutamate receptor, metabotropic 1 (Grm1), transcript variant 2
Database Link:	Entrez Gene 24414 Rat
Synonyms:	Metabotropic glutamate receptor 1, Gprc1a, Mglur1

Product images:





Western blot detection of mGluR1
20 ug of human brain lysate (lane 2)
with GluR1 polyclonal at 1:1000
dilution followed by AP-conjugated
secondary at 1:5000 dilution.
MW marker in Lane 1.



Western blot detection of mGluR1
20 ug of human brain lysate (lane 2)
with GluR1 polyclonal at 1:1000
dilution followed by AP-conjugated
secondary at 1:5000 dilution.
MW marker in Lane 1.