


Product datasheet for RA20064-100

TAC1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Mammalian, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens tachykinin precursor 1 (TAC1), transcript variant beta
Database Link:	Entrez Gene 24806 Rat
Synonyms:	Substance-P
Protein Families:	Druggable Genome, Secreted Protein


Product images:



Detection of Substance P in rat spinal cord dorsal horn (25 um sections). Courtesy of Chris Langdale, Urogenix.



[View online »](#)



Detection of Substance P in rat spinal cord dorsal horn (25 um sections). Courtesy of Chris Langdale, Urogenix.