


Product datasheet for RA24287-100

NOS1 (1419-1433) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Reactivity:	Feline, Guinea Pig, Human, Mouse, Primate, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	C-terminal synthetic peptide sequence corresponding to amino acids (1419-1433) of human nNos coupled to KLH.
Formulation:	State: Serum State: Lyophilized. 100 ul with <0.09% sodium azide as a preservative.
Gene Name:	nitric oxide synthase 1
Database Link:	Entrez Gene 4842 Human
Synonyms:	bNOS, nNOS, NOS type I, Neuronal NOS, Constitutive NOS, NC-NOS


Product images:




nNOS staining of rat brain.



[View online »](#)



nNOS staining of rat brain



nNOS staining of rat brain.