

## Product datasheet for **RA24288-100**

### 5 HT 2A (HTR2A) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens 5-hydroxytryptamine (serotonin) receptor 2A, G protein-coupled (HTR2A), transcript variant 1
Database Link:	<a href="#">Entrez Gene 3356 Human</a>
Synonyms:	5-hydroxytryptamine receptor 2A, 5-HT-2A, 5-HT-2
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Calcium signaling pathway, Gap junction, Neuroactive ligand-receptor interaction



[View online »](#)