


Product datasheet for RA24330-100

Serotonin transporter (SLC6A4) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens solute carrier family 6 member 4 (SLC6A4)
Database Link:	Entrez Gene 25553 Rat
Synonyms:	5HT transporter, 5HTT, HTT, Sodium-dependent serotonin transporter
Protein Families:	Druggable Genome, Transmembrane

Product images:



[View online »](#)

