

Product datasheet for RA24505-100

HTR2C Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens 5-hydroxytryptamine (serotonin) receptor 2C, G protein-coupled (HTR2C), transcript variant 2
Database Link:	Entrez Gene 15560 Mouse Entrez Gene 25187 Rat Entrez Gene 3358 Human
Synonyms:	HTR1C, 5HT-1C, 5-hydroxytryptamine receptor 2C, 5-HTR2C, 5-HT-2C, 5-HT2C
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Calcium signaling pathway, Gap junction, Neuroactive ligand-receptor interaction



[View online »](#)