

Product datasheet for **RA25004-100**

Somatostatin Receptor 1 (SSTR1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IP, WB
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens somatostatin receptor 1 (SSTR1)
Database Link:	Entrez Gene 20605 Mouse Entrez Gene 25033 Rat Entrez Gene 6751 Human
Synonyms:	Somatostatin receptor type 1, SS-1-R, SS1-R, SS1R, SRIF-2
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction



[View online »](#)