


## Product datasheet for RA25024-100

### EGR2 (200-300) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens early growth response 2 (EGR2), transcript variant 1
Database Link:	<a href="#">Entrez Gene 1959 Human</a>
Synonyms:	Early growth response protein 2, EGR-2, AT591, KROX20, Krox 20
Protein Families:	Druggable Genome, Transcription Factors


#### Product images:



Detection of EGR2 in human fetal lung tissue.



[View online »](#)



Detection of EGR2 in human fetal lung tissue.