


Product datasheet for RA25044-100

FOXD3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens forkhead box D3 (FOXD3)
Database Link:	Entrez Gene 15221 Mouse Entrez Gene 27022 Human
Synonyms:	Forkhead box protein D3


Product images:



Western blot analysis of FoxD3 (R233) antibody in extracts from HUVEC cells.



[View online »](#)



Western blot analysis of FoxD3 (R233) antibody in extracts from HUVEC cells.