


Product datasheet for RA25045-100

Oct4 (POU5F1) (100-200) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human, Mouse, Primate
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	POU class 5 homeobox 1
Database Link:	Entrez Gene 5460 Human
Synonyms:	POU5F1, OCT3, OCT4, OTF3, POU domain, class 5, transcription factor 1, Octamer-binding transcription factor 3, Oct-3, Oct-4


Product images:



Detection of OCT4 in mouse brain lysate.



[View online »](#)



Detection of OCT4 in mouse brain lysate.