


## Product datasheet for RA25066-100

### HIF-1 alpha (HIF1A) (400-550) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human, Primate
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	hypoxia inducible factor 1 alpha subunit
Database Link:	<a href="#">Entrez Gene 3091 Human</a>
Synonyms:	HIF-1 alpha, Hypoxia-inducible factor 1 alpha, ARNT-interacting protein, Member of PAS protein 1, Basic-helix-loop-helix-PAS protein MOP1, BHLHE78


#### Product images:



Detection of HIF-1 alpha. Lane 1: COS7 CoCl<sub>2</sub> treated cells Lane 2: COS7 untreated cells



[View online »](#)



Detection of HIF-1 alpha. Lane 1: COS7 CoCl<sub>2</sub> treated cells Lane 2: COS7 untreated cells