


## Product datasheet for RA25068-100

### LDL Receptor (LDLR) (500-550) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human, Mouse, Primate
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens low density lipoprotein receptor (LDLR), transcript variant 1
Database Link:	<a href="#">Entrez Gene 16835 Mouse</a> <a href="#">Entrez Gene 3949 Human</a>
Synonyms:	LDLR, Low-density lipoprotein receptor, LDL receptor
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Transmembrane
Protein Pathways:	Endocytosis


#### Product images:



Detection of LDL R by Western Blot



[View online »](#)



Detection of LDL R by Western Blot