


## Product datasheet for RA25069-100

### Nbn (350-400) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Mus musculus nibrin (Nbn)
Database Link:	<a href="#">Entrez Gene 4683 Human</a> <a href="#">Entrez Gene 27354 Mouse</a>
Synonyms:	NBN, NBS, NBS1, P95


#### Product images:




NBS1 staining of paraffin-embedded mouse bone marrow.




[View online »](#)



NBS1 Antibody-Western Blot on NIH/3T3 (Lane 4) and HeLa whole cell extract (Lane 5).



NBS1 staining of paraffin-embedded mouse bone marrow.



NBS1 Antibody-Western Blot on NIH/3T3 (Lane 4) and HeLa whole cell extract (Lane 5).