


Product datasheet for RA25070-100

TLR9 (1050-1100) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human, Primate, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens toll like receptor 9 (TLR9)
Database Link:	Entrez Gene 338457 Rat Entrez Gene 54106 Human
Synonyms:	Toll-like receptor 9, UNQ5798/PRO19605
Protein Families:	Druggable Genome, Transmembrane
Protein Pathways:	Toll-like receptor signaling pathway


Product images:



TLR9 Western Blot on Jurkat whole cell lysate.
Protocol on data-sheet.



[View online »](#)



TLR9 Western Blot on Jurkat whole cell lysate.
Protocol on data-sheet.