


## Product datasheet for RA25075-100

### Hif3a (N-term) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Mus musculus hypoxia inducible factor 3, alpha subunit (Hif3a), transcript variant 2
Database Link:	<a href="#">Entrez Gene 53417 Mouse</a>
Synonyms:	HIF-3 alpha, Hypoxia-inducible factor 3 alpha, BHLHE17, MOP7, PASD7, Member of PAS protein 7, IPAS


#### Product images:



Western blot analysis of HIF-3 alpha in Jurkat cell.



[View online »](#)



Western blot analysis of HIF-3 alpha in Jurkat cell.