


## Product datasheet for RA25081-100

### Olig2 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, IP, WB
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Mus musculus oligodendrocyte transcription factor 2 (Olig2)
Database Link:	<a href="#">Entrez Gene 10215 Human</a> <a href="#">Entrez Gene 304103 Rat</a> <a href="#">Entrez Gene 50913 Mouse</a>
Synonyms:	OLIG2, BHLHB1, PRKCBP2, RACK17, Oligodendrocyte transcription factor 2, Oligo2, bHLHB1

#### Product images:



[View online »](#)