

Product datasheet for RA25083-100

Bmi1 (COMMD3-BMI1) (1-100) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, WB
Reactivity:	Bovine, Canine, Chicken, Feline, Human, Mouse, Primate, Rabbit, Rat, Zebrafish
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens COMMD3-BMI1 readthrough (COMMD3-BMI1)
Database Link:	Entrez Gene 100532731 Human
Synonyms:	BMI1, PCGF4


Product images:



Bmi1 staining in HeLa cells.




[View online »](#)



Detection of Bmi1 protein in skeletal muscle.



Bmi1 staining in Hela cells.



Detection of Bmi1 protein in skeletal muscle.