


## Product datasheet for **RA25086-100**

### **GAP43 Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	FC, IHC, WB
Reactivity:	Chicken, Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens growth associated protein 43 (GAP43), transcript variant 2
Database Link:	<a href="#">Entrez Gene 14432 Mouse</a> <a href="#">Entrez Gene 29423 Rat</a> <a href="#">Entrez Gene 2596 Human</a>
Synonyms:	Neuromodulin


#### **Product images:**



GAP43 staining of region CA1-CA4 of the hippocampus. Note intense staining of cell processes.



[View online »](#)



GAP43 staining of region CA1-CA4 of the hippocampus. Note intense staining of cell processes.