


## Product datasheet for RA25095-100

### JARID2 (1-100) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ChIP, IF, WB
Reactivity:	Human, Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	jumonji and AT-rich interaction domain containing 2
Database Link:	<a href="#">Entrez Gene 3720 Human</a>
Synonyms:	JMJ


#### Product images:




JARID2 staining of HeLa cells using a DyLight 488 conjugated secondary antibody.




[View online »](#)




JARID2 staining of HeLa cells using a DyLight 488 conjugated secondary antibody




Western blot detection of JARID2 protein in E16 mouse brain lysates.



Western blot detection of JARID2 protein in E16 mouse brain lysates.



JARID2 staining of HeLa cells using a DyLight 488 conjugated secondary antibody.



Western blot detection of JARID2 protein in E16 mouse brain lysates.