


Product datasheet for RA26000-100

Angiopoietin 1 (ANGPT1) (C-term) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Reactivity:	Canine, Human, Rabbit, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens angiopoietin 1 (ANGPT1), transcript variant 1
Database Link:	Entrez Gene 89807 Rat Entrez Gene 284 Human
Synonyms:	ANG-1, ANG1, ANGPT1, KIAA0003
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Secreted Protein


Product images:



Top: ANG-1 staining of paraffin-embedded rat lung tissue (smooth muscle of blood vessel) using.
Bottom: Detection of ANG-1 from rat cardiac tissue lysate in Western blot assay.



[View online »](#)



Top: ANG-1 staining of paraffin-embedded rat lung tissue (smooth muscle of blood vessel) using. Bottom: Detection of ANG-1 from rat cardiac tissue lysate in Western blot assay.