


Product datasheet for RA26001-100

Mitofusin 2 (MFN2) (C-term) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Reactivity:	Human, Mouse, Rabbit, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens mitofusin 2 (MFN2), transcript variant 2
Database Link:	Entrez Gene 64476 Rat Entrez Gene 170731 Mouse Entrez Gene 9927 Human
Synonyms:	MFN2, CPRP1, KIAA0214
Protein Families:	Transmembrane


Product images:



Top: MFN2 staining of paraffin-embedded rat cardiac muscle tissue. Bottom: Detection of MFN2 from rat cardiac muscle tissue lysate by Western blot.



[View online »](#)



Top: MFN2 staining of paraffin-embedded rat cardiac muscle tissue. Bottom: Detection of MFN2 from rat cardiac muscle tissue lysate by Western blot.