



Product datasheet for RA26003-100

PLGF (PGF) (N-term) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens placental growth factor (PGF), transcript variant 1
Database Link:	Entrez Gene 18654 Mouse Entrez Gene 94203 Rat Entrez Gene 5228 Human
Synonyms:	PGFL, PLGF, PIGF
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Bladder cancer, Focal adhesion, mTOR signaling pathway, Pancreatic cancer, Pathways in cancer, Renal cell carcinoma



Product images:



PIGF staining of paraffin-embedded human placental tissue. Bottom: Detection of PLGF from rat lung tissue lysate in Western blot assay



[View online »](#)



PLGF staining of paraffin-embedded human placental tissue. Bottom: Detection of PLGF from rat lung tissue lysate in Western blot assay