



Product datasheet for RA26005-100

VEGFD (C-term) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Reactivity:	Human, Mouse, Rabbit, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens c-fos induced growth factor (vascular endothelial growth factor D) (FIGF)
Database Link:	Entrez Gene 2277 Human
Synonyms:	FIGF, VEGF-D, FIGF
Protein Families:	Druggable Genome, Secreted Protein
Protein Pathways:	Bladder cancer, Cytokine-cytokine receptor interaction, Focal adhesion, mTOR signaling pathway, Pancreatic cancer, Pathways in cancer, Renal cell carcinoma


Product images:



Immunohistochemical staining of paraffin-embedded human breast cancer tissue, using Anti-VEGF-D. Bottom: Detection of VEGF-D from rat cardiac muscle tissue lysate in Western blot assay, using Anti-VEGF-D.



[View online »](#)



Immunohistochemical staining of paraffin-embedded human breast cancer tissue, using Anti-VEGF-D. Bottom: Detection of VEGF-D from rat cardiac muscle tissue lysate in Western blot assay, using Anti-VEGF-D.