


Product datasheet for RA30000-50

5HT1E Receptor (HTR1E) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Reactivity:	Human, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens 5-hydroxytryptamine receptor 1E (HTR1E)
Database Link:	Entrez Gene 3354 Human
Synonyms:	5-hydroxytryptamine receptor 1E, 5-HT-1E, 5-HT1E
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction


Product images:



5HT-1E staining of paraffin embedded testis.



[View online »](#)



5HT-1E staining of paraffin embedded testis.