


Product datasheet for RA30001-50

5HT5A receptor (HTR5A) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Reactivity:	Human, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens 5-hydroxytryptamine receptor 5A (HTR5A)
Database Link:	Entrez Gene 25689 Rat Entrez Gene 3361 Human
Synonyms:	5-hydroxytryptamine receptor 5A, 5-HT5A, 5-HT-5A, 5-HT-5
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Calcium signaling pathway, Neuroactive ligand-receptor interaction


Product images:



5HT5A Staining of paraffin embedded human cerebellum.



[View online »](#)



5HT5A Staining of paraffin embedded human cerebellum.