


Product datasheet for RA30002-50

5 HT1F (HTR1F) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens 5-hydroxytryptamine receptor 1F (HTR1F)
Database Link:	Entrez Gene 3355 Human
Synonyms:	5-hydroxytryptamine receptor 1F, HTR1EL, 5-HT-1F
Protein Families:	Druggable Genome, GPCR, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction


Product images:



5HT5A Staining of paraffin embedded human cerebellum.



[View online »](#)



5HT5A Staining of paraffin embedded human cerebellum.