


Product datasheet for RA30017-50

GFR alpha 3 (GFRA3) (347-360) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC, WB
Reactivity:	Mouse
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens GDNF family receptor alpha 3 (GFRA3)
Database Link:	Entrez Gene 14587 Mouse
Synonyms:	GFR-alpha-3, GDNF receptor alpha-3
Protein Families:	Druggable Genome


Product images:



GFRa-3 staining of Mouse Adrenal Tissue (FFPE)



[View online »](#)



GFRA-3 staining of Mouse Adrenal Tissue (FFPE)