


Product datasheet for **RA30023-50**

PACRG (204-215) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, IHC, WB
Reactivity:	Human, Rat, Zebrafish
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	PARK2 coregulated
Database Link:	Entrez Gene 135138 Human
Synonyms:	Parkin coregulated gene protein, GLUP, PARK2 coregulated gene protein, PARK2CRG


Product images:



PACRG Staining of human colon: Formalin-Fixed Paraffin-Embedded (FFPE).



[View online »](#)



PACRG Staining of human colon: Formalin-Fixed Paraffin-Embedded (FFPE).