


## Product datasheet for RA30035-50

### BMP4 (20-34) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Reactivity:	Bovine, Chicken, Frog, Human, Mouse, Porcine, Rat, Zebrafish
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	bone morphogenetic protein 4
Database Link:	<a href="#">Entrez Gene 652 Human</a>
Synonyms:	BMP2B, BMP-2B, BMP-4, DVR4, Bone morphogenetic protein 4


#### Product images:



BMP4 staining of Human Lung Tissue, respiratory epithelium: Formalin-Fixed Paraffin-Embedded (FFPE).



[View online »](#)



BMP4 staining of Human Lung Tissue, respiratory epithelium: Formalin-Fixed Paraffin-Embedded (FFPE).