


Product datasheet for RA30047-50

GABRP Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	gamma-aminobutyric acid type A receptor pi subunit
Database Link:	Entrez Gene 2568 Human
Synonyms:	GABA A receptor subunit pi


Product images:



GABRP expression in the Brain (Cerebellum):
Formalin-Fixed Paraffin-Embedded (FFPE).



[View online »](#)



GABRP expression in the Brain (Cerebellum):
Formalin-Fixed Paraffin-Embedded (FFPE).