


Product datasheet for RA30049-50

GPRC6A (N-term) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Reactivity:	Human
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens G protein-coupled receptor class C group 6 member A (GPRC6A), transcript variant 1
Database Link:	Entrez Gene 222545 Human
Synonyms:	G-protein coupled receptor family C group 6 member A
Protein Families:	Druggable Genome, Transmembrane


Product images:



GPRC6A of Human Skeletal Muscle (paraffin-embedded).



[View online »](#)



GPRC6A of Human Skeletal Muscle (paraffin-embedded).