


Product datasheet for RA30054-50

TLR10 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens toll like receptor 10 (TLR10), transcript variant 2
Database Link:	Entrez Gene 81793 Human
Synonyms:	Toll-like receptor 10, UNQ315/PRO358
Protein Families:	Druggable Genome, Transmembrane


Product images:



TLR2 staining of Human Lung, Alveolar Macrophages (formalin-fixed, paraffin-embedded) at 5 ug/ml followed by biotinylated goat anti-rabbit IgG secondary antibody, alkaline phosphatase-streptavidin and chromogen. Protocol on data-sheet.



[View online »](#)



TLR2 staining of Human Lung, Alveolar Macrophages (formalin-fixed, paraffin-embedded) at 5 ug/ml followed by biotinylated goat anti-rabbit IgG secondary antibody, alkaline phosphatase-streptavidin and chromogen. Protocol on data-sheet.