

Product datasheet for **RC200008A1V**

Human **TEX264 (NM_015926) AAV Particle**

Product data:

Product Type: AAV Particles
Product Name: Human TEX264 (NM_015926) AAV Particle
Tag: Myc-DDK
Symbol: TEX264
Synonyms: ZSIG11
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC200008 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGTCGGACCTGCTACTACTGGGCCTGATTGGGGCCTGACTCTCTTACTGCTGCTGACGCTGCTGGCCT
 TTGCCGGTACTCAGGGCTACTGGCTGGGGTGAAGTGAGTGCTGGGTCACCCCCATCCGCAACGTCAC
 TGTGGCCTACAAGTCCACATGGGGCTCTATGGTGAGACTGGGCGCTTTTCACTGAGAGCTGCAGCATC
 TCTCCAAGCTCCGCTCCATCGCTGTCTACTATGACAACCCCCACATGGTGCCCTGATAAGTGCCGAT
 GTGCCGTGGGCAGCATCCTGAGTGAAGGTGAGGAATCGCCCTCCCTGAGCTCATCGACCTTACCAGAA
 ATTTGGCTTCAAGGTGTTCTCCTTCCCGCACCCAGCCATGTGGTGACAGCCACCTTCCCTACACCACC
 ATTCTGTCCATCTGGCTGGCTACCCGCGTGTCCATCCTGCCTTGGACACCTACATCAAGGAGCGGAAGC
 TGTGTGCCTATCCTCGGCTGGAGATCTACCAGGAAGACCAGATCCATTTTCATGTGCCACTGGCACGGCA
 GGGAGACTTCTATGTGCCTGAGATGAAGGAGACAGAGTGGAAATGGCGGGGGCTTGTGGAGGCCATTGAC
 ACCCAGGTGGATGGCACAGGAGCTGACACAATGAGTGACACGAGTTCTGTAAGCTTGAAGTGAGCCCTG
 GCAGCCGGGAGACTTCACTGCCACACTGTACCTGGGCGAGCAGCCGTGGCTGGGATGACGGTGACAC
 CCGCAGCGAGCACAGCTACAGCGAGTCAGGTGCCAGCGCTCCTCTTTGAGGAGCTGGACTTGGAGGGC
 GAGGGGCCCTTAGGGGAGTCACGGCTGGACCCTGGGACTGAGCCCTGGGGACTACCAAGTGCTCTGGG
 AGCCCACTGCCCTGAGAAGGGCAAGGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	<p>>RC200008 protein sequence Red=Cloning site Green=Tags(s)</p> <p>MSDLLLLGLIGGLTLLLLLTLAFAGYSGLLAGVEVSAGSPPIRNVTVA YKFHMGLYGETGRLFTESCSI SPKLRSIAVYYDNPHMVPPDKCRCAVGSILSEGEESPELIDLYQKFGKVF SFPAPSHVVTATFPYTT ILSIWLATRRVHPALD TYIKERKLCAYPRLEIYQEDQIHFMCLARQGDFYVP EMKETEWKWRGLVEAID TQVDGTGADTMSDTSSVSLEVSPG SRETSAA TLSPGASSRGWDDGDTRSEHSYSESGASGSSFEELDLEG EGPLGESRLDPGTEPLGTTKWLWEPTAPEKGKE</p> <p>TRTRPLEQKLI SEEDLAANDILDYKDDDDKV</p>
Species:	Human
Serotype:	AAV-2
ACCN:	NM_015926
ORF Size:	939 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_015926.3
RefSeq Size:	1403 bp
RefSeq ORF:	942 bp
Locus ID:	51368
UniProt ID:	Q9Y6I9
Cytogenetics:	3p21.2
MW:	34.2 kDa