

## Product datasheet for **RC200009L1V**

### PPME1 (NM\_016147) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PPME1 (NM_016147) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PPME1
Synonyms:	ABDH19; PME-1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_016147
ORF Size:	1158 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200009).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_016147.1</a>
RefSeq Size:	2484 bp
RefSeq ORF:	1161 bp
Locus ID:	51400
UniProt ID:	<a href="#">Q9Y570</a>
Cytogenetics:	11q13.4
Protein Families:	Druggable Genome
MW:	42.1 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein phosphatase methylesterase localized to the nucleus. The encoded protein acts on the protein phosphatase-2A catalytic subunit and supports the ERK pathway through dephosphorylation of regulatory proteins. It plays a role in malignant glioma progression. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2012]