

Product datasheet for **RC200020A1V**

Human TPPP3 (NM_016140) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human TPPP3 (NM_016140) AAV Particle
Tag: Myc-DDK
Symbol: TPPP3
Synonyms: CGI-38; p20; p25gamma; TPPP/p20
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC200020 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCAGCGAGCACAGACATGGCTGGGCTGGAGGAGAGCTTCCGCAAGTTGCCATCCATGGTGACCCCA
AGGCCAGTGGGCAAGAGATGAATGGCAAGAACTGGGCCAAGCTGTGCAAGGACTGCAAGGTGGCTGACGG
AAAGTCCGTGACAGGGACCGATGTGGACATCGTCTTCTCCAAAGTCAAGGGGAAGTCTGCTCGGGTCATC
AACTATGAGGAGTTCAAGAAGGCCCTGGAAGAGCTGGCGACCAAGAGATTCAAGGGGAAGAGCAAGGAGG
AGGCCTTCGATGCCATCTGCCAGCTGGTGGCAGGCAAGAGCCAGCCAATGTGGGCGTCACTAAAGCAAA
AACAGGGGGTGTCTGTAGACCGGCTGACGGACACCAGCAGATACACGGGGCTCCCAAGGAGCGCTTCGAT
GAGAGCGGCAAGGGCAAGGGCATTGCGGGACGGCAGGACATCCTGGACGACAGTGGCTACGTGAGCGCCT
ACAAGAATGCAGGCACCTACGATGCCAAGGTGAAGAAG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC200020 protein sequence
Red=Cloning site Green=Tags(s)

MAASTDMAGLEESFRKFAIHGDPKASGQEMNGKNWAKLCKDCKVADGKSVTGTVDIVFSKVKGKSARVI
NYEEFKKALEELATKRFKGSKEEAFDAICQLVAGKEPANVGVTAKTGGAVDRLTDSRYTGSHKERFD
ESGKGKGIAGRQDILDDSGYVSAYKNAGTYDAKVKK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_016140
ORF Size:	528 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_016140.2</u>
RefSeq Size:	1159 bp
RefSeq ORF:	531 bp
Locus ID:	51673
UniProt ID:	<u>Q9BW30</u>
Cytogenetics:	16q22.1
MW:	19 kDa