

## Product datasheet for **RC200033L1V**

### Mitochondrial ribosomal protein L11 (MRPL11) (NM\_016050) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Mitochondrial ribosomal protein L11 (MRPL11) (NM_016050) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Mitochondrial ribosomal protein L11
Synonyms:	CGI-113; L11MT; MRP-L11
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_016050
ORF Size:	576 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200033).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_016050.2</a>
RefSeq Size:	1622 bp
RefSeq ORF:	579 bp
Locus ID:	65003
UniProt ID:	<a href="#">Q9Y3B7</a>
Cytogenetics:	11q13.2
Domains:	Ribosomal_L11


[View online »](#)

**MW:** 20.7 kDa

**Gene Summary:** This nuclear gene encodes a 39S subunit component of the mitochondrial ribosome. Alternative splicing results in multiple transcript variants. Pseudogenes for this gene are found on chromosomes 5 and 12. [provided by RefSeq, May 2014]