

Product datasheet for **RC200042L3V**

Bcl rambo (BCL2L13) (NM_015367) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Bcl rambo (BCL2L13) (NM_015367) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BCL2L13
Synonyms:	BCL-RAMBO; Bcl2-L-13; MIL1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_015367
ORF Size:	1455 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200042).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015367.2
RefSeq Size:	5132 bp
RefSeq ORF:	1458 bp
Locus ID:	23786
UniProt ID:	Q9BXK5
Cytogenetics:	22q11.21
Domains:	Bcl-2
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

MW: 52.7 kDa

Gene Summary: This gene encodes a mitochondrially-localized protein with conserved B-cell lymphoma 2 homology motifs. Overexpression of the encoded protein results in apoptosis. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Jul 2012]