

Product datasheet for **RC200053L3V**

POLDIP2 (NM_015584) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	POLDIP2 (NM_015584) Human Tagged ORF Clone Lentiviral Particle
Symbol:	POLDIP2
Synonyms:	p38; PDIP38; POLD4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_015584
ORF Size:	1104 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200053).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015584.3
RefSeq Size:	2753 bp
RefSeq ORF:	1107 bp
Locus ID:	26073
UniProt ID:	Q9Y2S7
Cytogenetics:	17q11.2
Domains:	DUF525
MW:	42 kDa


[View online »](#)

Gene Summary:

This gene encodes a protein that interacts with the DNA polymerase delta p50 subunit, as well as with proliferating cell nuclear antigen. The encoded protein maybe play a role in the ability of the replication fork to bypass DNA lesions. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2014]