

Product datasheet for **RC200076**

DBNDD2 (NM_001048221) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: DBNDD2 (NM_001048221) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: DBNDD2
Synonyms: C20orf35; CK1BP; HSMNP1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC200076 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACCCAAATCCTCGGGCCGCCCTGGAGCGCCAGCAGCTCCGCCTTCGGGAGCGGCAAAAATTCTTCG
AGGACATTTTACAGCCAGAGACAGAGTTTGTCTTCTCTGTCCCATCTGCATCTCGAGTCGAGAGACC
CCCCATAGGTAGTATCTCATCCATGGAAGTGAATGTGGACACACTGGAGCAAGTAGAACTTATTGACCTT
GGGACCCGGATGCAGCAGATGTGTTCTTGCCTTGCGAAGATCCTCCACCAACCCCCAGTCGTCTGGG
TGGACAACCATTTGGAGGAGCTGAGCCTGCCGGTGCCTACATCAGACAGGACCACATCTAGGACCTCCTC
CTCCTCCTCCTCCGACTCCTCCACCAACCTGCATAGCCCAAATCCAAGTGATGATGGAGCAGATACGCC
TTGGCACAGTCGGATGAAGAGGAGAAAGGGGTGATGGAGGGGCAGAGCCTGGAGCTGCAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC200076 protein sequence
Red=Cloning site Green=Tags(s)
MDPNPRAALERQQLRLRERQKFFEDILQPETEFVPLSHLHLESQRPPIGSISSMEVNVDTLEQVELIDL
GDPDAADVFLPCEDPPPTPQSSGVDNHLEELSLPVPTSDRTTSRTSSSSSSDSTNLHSPNPSDDGADTP
LAQSDEEEERGDGGAEPGACS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6384_e10.zip



[View online >](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001048221

ORF Size: 483 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001048221.3](#)


RefSeq Size: 1165 bp

RefSeq ORF: 486 bp

Locus ID: 55861

UniProt ID: [Q9BQY9](#)
Cytogenetics: 20q13.12
MW: 17.5 kDa
Gene Summary: May modulate the activity of casein kinase-1. Inhibits CSNK1D autophosphorylation (in vitro). [UniProtKB/Swiss-Prot Function]


Product images:



Circular map for RC200076



Coomassie blue staining of purified DBNDD2 protein (Cat# [TP300076]). The protein was produced from HEK293T cells transfected with DBNDD2 cDNA clone (Cat# RC200076) using MegaTran 2.0 (Cat# [TT210002]).



Western blot validation of overexpression lysate (Cat# [LY420788]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with [RC214187] using transfection reagent MegaTran 2.0 (Cat# [TT210002]).