

## Product datasheet for **RC200079L2V**

### SCAND1 (NM\_033630) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SCAND1 (NM_033630) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SCAND1
Synonyms:	RAZ1; SDP1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_033630
ORF Size:	537 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200079).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_033630.1</a> , <a href="#">NP_361012.1</a>
RefSeq Size:	919 bp
RefSeq ORF:	540 bp
Locus ID:	51282
UniProt ID:	<a href="#">P57086</a>
Cytogenetics:	20q11.23
Domains:	LER
Protein Families:	Transcription Factors



[View online »](#)

**MW:** 19.1 kDa

**Gene Summary:** This gene encodes a SCAN box domain-containing protein. The SCAN domain is a highly conserved, leucine-rich motif of approximately 60 aa originally found within a subfamily of zinc finger proteins. This gene belongs to a family of genes that encode an isolated SCAN domain, but no zinc finger motif. This protein binds to and may regulate the function of the transcription factor myeloid zinc finger 1B. Alternate splicing results in multiple transcript variants.[provided by RefSeq, Jan 2011]