

Product datasheet for **RC200088A1V**

Human FIBP (NM_004214) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human FIBP (NM_004214) AAV Particle
Tag:	Myc-DDK
Symbol:	FIBP
Synonyms:	FGFIBP; FIBP-1; TROFAS
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200088 representing NM_004214 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGACCACTGAGCTGGACATCTTCGTGGGGAACACGACCCTTATCGACGAGGACGTGTATCGCCTCTGGC
TCGATGGTTACTCGGTGACCGACGCGGTGGCCCTGCGGGTGGCTCGGGAATCCTGGAGCAGACTGGCGC
CACGGCAGCGGTGCTGCAGAGCGACACCATGGACCATTACCGCACCTCCACATGCTCGAGCGGTGCTG
CATGCGCCGCCAAGCTACTGCACCACTCATCTTCAGATTCCGCCCTCCCGGAGGCACTACTCATCG
AGAGGTACTATGCCTTTGATGAGGCCTTTGTTCCGGAGGTGCTGGCAAGAAGCTGTCCAAGGCACCAA
GAAAGACCTGGATGACATCAGCACCAAAACAGGCATCACCTCAAGAGCTGCCGGAGACAGTTTGACAAC
TTTAAACGGGTCTTCAAGGTGGTAGAGGAAATGCGGGGCTCCCTGGTGGACAATATTCAGCAACACTTCC
TCCTCTCTGACCGGTTGGCCAGGGACTATGCAGCCATCGTCTTCTTTGCTAACAACCGCTTTGAGACAGG
GAAGAAAAAAGTGCAGTATCTGAGCTTCGGTGACTTTGCCTTCTGCGCTGAGCTCATGATCCAAAAGTGG
ACCTTTGGAGCCGTCGACTCACAGATGGATGACATGGACATGGACTTAGACAAGGAATTTCTCCAGGACT
TGAAGGAGCTCAAGGTGCTAGTGGCTGACAAGGACCTTCTGGACCTGCACAAGAGCCTGGTGTGCACTGC
TCTCCGGGGAAGCTGGCGCTCTCTGAGATGGAAGCAACTCAAGAACCTGTCCCGGGGGCTGGTG
AACGTGGCCGCCAAGCTGACCCACAATAAGATGTGAGAGACCTGTTTGTGGACCTCGTGGAGAAGTTTG
TGGAACCTGCGCTCCGACCACTGGCCACTCAGCGACGTGCGGTTCTTCCTGAATCAGTATTCAGCGTC
TGTCCTCTCCCTCGATGGCTTCGACACCCAGGCCCTCTGGGACCGCTACATGGGACCCCTCCGCGGCTGC
CTCCTGCGCCTGTATCATGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200088 representing NM_004214 Red=Cloning site Green=Tags(s) MTSELDIFVGNTTLIDEDVYRLWLDGYSVTDAVALRVRSGILEQTGATAAVLQSDTMDHYRTFHMLERLL HAPPKLLHQLIFQIPPSRQALLIERYYAFDEAFVREVLGKKLSKGTKKDLDDISTKTGITLKSCRRQFDN FKRVFKVVEEMRGSLVDNIQQHFLLSDRLARDYAAIVFFANNRFETGKKKLQYLSFGDFAFCAELMIQNW TLGAVDSQMDDMDMDLDKEFLQDLKELKVLVADKDLLDLHKS LVCTALRGKLGVFSEMEANFKNLSRGLV NYAAKLTHNKDVRDLFVDLVEKFVEPCRS DHWPLSDVRFFLNQYSASVHSLDGFRHQALWDRYMGTLRGC LLRLYHD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_004214
ORF Size:	1071 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_004214.4</u>
RefSeq Size:	1429 bp
RefSeq ORF:	1074 bp
Locus ID:	9158
UniProt ID:	<u>O43427</u>
Cytogenetics:	11q13.1
MW:	41 kDa