

Product datasheet for **RC200090L4V**

RPAP1 (NM_015540) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | RPAP1 (NM_015540) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | RPAP1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_015540 |
| ORF Size: | 4179 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC200090). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_015540.2 |
| RefSeq Size: | 4698 bp |
| RefSeq ORF: | 4182 bp |
| Locus ID: | 26015 |
| UniProt ID: | Q9BWH6 |
| Cytogenetics: | 15q15.1 |
| MW: | 152.6 kDa |
| Gene Summary: | This protein forms part of the RNA polymerase II (RNAPII) enzyme complex and may recruit RNAPII to chromatin through its interaction with acetylated histones. [provided by RefSeq, Jul 2012] |


[View online »](#)