

## Product datasheet for **RC200101L3V**

### SLC35C1 (NM\_018389) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC35C1 (NM_018389) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC35C1
Synonyms:	CDG2C; FUCT1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018389
ORF Size:	1092 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200101).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018389.3</a> , <a href="#">NP_060859.3</a>
RefSeq Size:	3531 bp
RefSeq ORF:	1095 bp
Locus ID:	55343
UniProt ID:	<a href="#">Q96A29</a>
Cytogenetics:	11p11.2
Protein Families:	Transmembrane
MW:	39.7 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a GDP-fucose transporter that is found in the Golgi apparatus. Mutations in this gene result in congenital disorder of glycosylation type IIc. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2009]