

Product datasheet for **RC200105A1V**

Human TRUB2 (NM_015679) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TRUB2 (NM_015679) AAV Particle
Tag:	Myc-DDK
Symbol:	TRUB2
Synonyms:	CLONE24922
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC200105 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGGTCTGCTGGCTTGTGCGGGTGCATGGGCTTTTCGCGGTCTATAAGCCCCGGGGCTAAATGGA
AGCACCTGCGGGATACAGTGGAGCTACAACCTCTGAAGGGTCTCAATGCCAGGAAGCCTCCCGCTCTAA
ACAGCGTGTTCGCTTCTTGCTGGGCCCCATGGAAGGCAGCGAAGAGAAGGAGCTGACCTCACAGCCACC
AGCGTACCCTCTTTCATCAACCATCCACTGGTATGTGGACCAGCATTGCGCCATCTCAAGGTTGGCGTGG
GACATCGGTTGGATGCCAGGCTTCTGGAGTACTTGTGCTCGGCGTGGGACATGGATGCAGGCTCCTCAC
CGATATGTACAATGCTCATCTTACCAAGGATTACACAGTGCCTGGCCTCCTGGGCAAAGCTACAGATGAC
TTCCGTGAGGACGGGAGGCTGGTAGAGAAGACAACCTATGACCACGTGACCAGAGAGAAGCTGGACCGCA
TTCTGGCCGTTATCCAAGGCTCCCATCAGAAGGCCCTGGTGTGTACTCCAACCTCGACCTGAAGACCCA
GGAGGCCTATGAGATGGCCGTGAGAGGCCTGATCCGGCCCATGAACAAGTCCCCGATGCTGATAACTGGC
ATCCGATGCCTCTACTTTGCACCTCCGGAATTCCTCTTAGAGGTGCAGTGCATGCATGAGACGCAGAAAG
AGCTGCGGAAGTTGGTTCATGAAATCGGCCTGGAATAAGACCACTGCTGTCTGCACCAAGTGGCGG
CACGCGCGACGGCTTCTTACGCTAGACAGTGCCTCCTGAGGACCCAGTGGGACCTAACCAACATCCAG
GATGCTATCCGGGCTGCTACCCCTCAGGTAGCTGCAGAGCTGGAGAAGAGCTTGAGCCCGGGGCTGGACA
CCAAGCAGCTCCCAAGTCCGGGATGGTCTGGGACTCCAGGGCCCGAGCTCTACCTTGGGGCTGGAGAG
GGGTGCGGGGACAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC200105 protein sequence Red=Cloning site Green=Tags(s) MGSAGLSRLHGLFAVYKPPGLKWKHLRDTVELQLLKGLNARKPPAPKQVRVFLGPMEGSEEKELTLTAT SVPSFINHPLVCGPAFAHLKVGVGHRDAQASGVLVLGVGHGCRLLTDMYNAHLTKDYTVRGLLGKATDD FREDGRLVEKTTYDHTREKLDRLAVIQGSHQKALVMYSNLDLKTQEAYEMAVRGLIRPMNKSPMLITG IRCLYFAPPEFLLEVQCMHETQKELRKLVEIGLELKTAVCTQVRRTRDGGFTLDSALLRTQWDLTNIQ DAIRAATPQVAAELEKSLSPGLDTKQLSPGWSWDSQGPSSTLGLERGAGQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_015679
ORF Size:	993 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_015679.1</u>
RefSeq Size:	1439 bp
RefSeq ORF:	996 bp
Locus ID:	26995
UniProt ID:	<u>O95900</u>
Cytogenetics:	9q34.11
MW:	36.7 kDa