

Product datasheet for **RC200109L4V**

DHRS4 (NM_021004) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DHRS4 (NM_021004) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DHRS4
Synonyms:	CR; NRDR; PHCR; PSCD; SCAD-SRL; SDR-SRL; SDR25C1; SDR25C2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021004
ORF Size:	693 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200109).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021004.2
RefSeq Size:	1315 bp
RefSeq ORF:	837 bp
Locus ID:	10901
UniProt ID:	Q9BTZ2
Cytogenetics:	14q11.2
Domains:	adh_short
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Metabolic pathways, Retinol metabolism

MW: 24.6 kDa

Gene Summary: Reduces all-trans-retinal and 9-cis retinal. Can also catalyze the oxidation of all-trans-retinol with NADP as co-factor, but with much lower efficiency. Reduces alkyl phenyl ketones and alpha-dicarbonyl compounds with aromatic rings, such as pyrimidine-4-aldehyde, 3-benzoylpyridine, 4-benzoylpyridine, menadione and 4-hexanoylpyridine. Has no activity towards aliphatic aldehydes and ketones (By similarity).[UniProtKB/Swiss-Prot Function]