

Product datasheet for RC200112L3

IMP3 (NM_018285) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: IMP3 (NM_018285) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: IMP3

BRMS2; C15orf12; MRPS4 Synonyms:

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

E. coli Selection: Chloramphenicol (34 ug/mL)


The ORF insert of this clone is exactly the same as(RC200112). **ORF Nucleotide**

Sequence:

Restriction Sites: Sgfl-RsrII

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_018285

ORF Size: 552 bp



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



IMP3 (NM_018285) Human Tagged Lenti ORF Clone - RC200112L3

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

15q24.2

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 018285.2</u>

RefSeq Size: 1255 bp RefSeq ORF: 555 bp

Locus ID: 55272

UniProt ID: Q9NV31

Domains: S4

Cytogenetics:

MW: 21.9 kDa


Gene Summary: This gene encodes the human homolog of the yeast Imp3 protein. The protein localizes to the

nucleoli and interacts with the U3 snoRNP complex. The protein contains an S4 domain.

[provided by RefSeq, Jul 2008]



Product images:



Circular map for RC200112L3