

Product datasheet for **RC200116L4V**

KLHDC2 (NM_014315) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	KLHDC2 (NM_014315) Human Tagged ORF Clone Lentiviral Particle
Symbol:	KLHDC2
Synonyms:	HCLP-1; HCLP1; LCP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014315
ORF Size:	1218 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC200116).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014315.2
RefSeq Size:	1808 bp
RefSeq ORF:	1221 bp
Locus ID:	23588
UniProt ID:	Q9Y2U9
Cytogenetics:	14q21.3
Domains:	Kelch
MW:	46.1 kDa



[View online »](#)

Gene Summary: Represses CREB3-mediated transcription by interfering with CREB3-DNA binding.
[UniProtKB/Swiss-Prot Function]